


Place cursor on number or click for description



- | | | |
|-------------------------------|----------------------------------|-----------------------------|
| 1 Donor line | 13 Return reservoir | 25 Inlet port |
| 2 Sample bag | 14 Return pump header | 26 Collection chamber |
| 3 AC/draw/return manifold | 15 Centrifuge pressure sensor | 27 LRS chamber |
| 4 Anticoagulant (AC) line | 16 Inlet line to centrifuge | 28 Plasma collect line |
| 5 Sterile barrier filter | 17 RBC line from centrifuge | 29 Plasma bag |
| 6 Cassette | 18 Platelet line from centrifuge | 30 Platelet collect line |
| 7 AC pump header | 19 Plasma line from centrifuge | 31 Platelet bag |
| 8 Draw/return pressure sensor | 20 Plasma pump header | 32 Platelet product sampler |
| 9 AC Line | 21 Cassette label | 33 Vent bag line |
| 10 Return line | 22 Collect pump header | 34 Vent bag |
| 11 Draw line | 23 Centrifuge loop | |
| 12 Inlet pump header | 24 Channel | |

Spectra Optia® Platelet, Plasma Set Catalog Numbers 10400, 10406

Spectra Optia[®] Platelet, Plasma Set

Part Descriptions

1. **Donor line:** provides access to and return to the donor.
2. **Sample bag:** used to collect blood from the donor for sampling purposes.
3. **AC/draw/return manifold:** consists of the access injection site and the connections for the anticoagulant line (4), the draw line (11), and the return line (10).
4. **Anticoagulant (AC) line (with orange spike):** carries anticoagulant from the AC bag to the cassette (6).
5. **Sterile barrier filter:** prevents bacteria from entering the system to maintain a functionally closed system.
6. **Cassette:** guides the flow of blood and products through the tubing set.
7. **AC pump header:** the tubing segment that fits into the AC pump.
8. **Draw/return pressure sensor:** monitors the draw and return pressures at the donor access site.
9. **AC line:** carries AC from the cassette (6) to the AC/draw/return manifold (3).
10. **Return line:** carries blood components back to the donor.
11. **Draw line:** carries anticoagulated whole blood into the tubing set.
12. **Inlet pump header:** the tubing segment that fits into the inlet pump.
13. **Return reservoir:** receives uncollected components for return to the donor. Also provides a return filter (200 micron) and air detection.
14. **Return pump header:** the tubing segment that fits into the return pump.
15. **Centrifuge pressure sensor:** detects high pressure in the centrifuge.
16. **Inlet line to centrifuge:** carries blood to the centrifuge.
17. **RBC line from centrifuge:** carries RBC (red blood cells) from the centrifuge for collection or return to the donor.
18. **Platelet line from centrifuge:** carries platelets from the centrifuge for collection or return to the donor.
19. **Plasma line from centrifuge:** carries plasma from the centrifuge for collection or return to the donor.
20. **Plasma pump header:** the tubing segment that fits into the plasma pump.
21. **Cassette label:** used by the Spectra Optia system RBC detector to identify a tubing set as capable of collecting platelet products.
22. **Collect pump header:** the tubing segment that fits into the collect pump.
23. **Centrifuge loop:** consists of lines that are loaded in the centrifuge. Includes:
 - Four-lumen tubing: used to carry fluid into and out of the channel.
 - Sleeves: used to reinforce the tubing.
 - Collars: used to secure the ends of the loop in the centrifuge.
 - Bearings: used to connect the loop to the centrifuge.
24. **Channel:** used during centrifugation to separate the donor's blood into cellular components.
25. **Inlet port:** routes incoming blood into the channel.
26. **Collection chamber:** routes separated blood products to the appropriate collection lines.
27. **LRS chamber:** provides leukoreduction of collected platelets.
28. **Plasma collect line:** carries the collected plasma to the plasma bag (29).
29. **Plasma bag:** 1 L bag that holds concurrently collected plasma.
30. **Platelet collect line:** carries the collected platelets to the platelet bags (31).
31. **Platelet bag:** the bag where the collected platelets are stored. The bag is made from PVC with a citrate plasticizer.
32. **Platelet product sampler:** used for product quality testing.
33. **Vent bag line:** carries displaced air to and from the vent bag (34).
34. **Vent bag:** holds displaced air from the system.